

ElectroForce® 3300 Series III



Site Preparation Guide

Table of Contents

Table of Contents	2
Ideal Setup	3
System Components.....	4
Instrument Measurements	5
Utility Requirements.....	7–8
Power	7
Gas	8
Hot/Cold Chamber.....	8
Site Preparation Checklist	9
TA Instrument Offices.....	10



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab



Customer

Ideal Setup



IDEAL PLACEMENT AND BENCH MEASUREMENTS

Select a location with adequate floor and ceiling space and a rigid laboratory bench that is level and is in a vibration-free environment. Bench must be rated to support several hundred pounds.



Distance from the wall:
0.15 m (0.5 ft) min.

Table width: 1.5 m (5 ft)


Table depth: 1.2 m (4 ft)

Floor space: 91.44 cm x
91.44 cm (3 ft x 3 ft)

System Components



MAIN SYSTEM COMPONENTS



- | | |
|------------------------|-------------------------|
| A. Test Instrument | E. Power Supply (Axial) |
| B. Computer Tower | F. Emergency Stop |
| C. Computer Monitor | |
| D. PCI Electronics Box | |

Instrument Measurements



3300 SERIES III – FLOOR STANDING FRAME



Height: 1854 mm (73 in)

Width: 579 mm (22.8 in)

Depth: 518 cm (20.4 in)

Weight: 279 kg (615 lbs)



3300 SERIES III – BENCHTOP

Height: 1232 mm (48.5 in)

Width: 579 mm (22.8 in)

Depth: 518 cm (20.4 in)

Weight: 172 kg (381 lbs)



Instrument Measurements



POWER SUPPLY FOR 3330

Height: 53.34 cm (21 in)

Width: 43.18 cm (17 in)

Depth: 40.64 cm (16 in)



POWER SUPPLY FOR 3310

Height: 30.48 cm (12 in)

Width: 34.29 cm (13.5 in)

Depth: 40.64 cm (16 in)



Utility Requirements



POWER

Item	Requirement
Instrument Power	<ul style="list-style-type: none"> • 200V, 50–60 Hz, 15A • 230V, 50–60 Hz, 13A • Neutral to Ground (NG) voltage max 0.5 volt • Safety ground per local regulation
PCI Box Power	<ul style="list-style-type: none"> • 104–120V, 50–60 Hz Hz, 1.2 A • 207–230V, 50–60 Hz Hz, 0.6 A
Torsion Power Supply	<ul style="list-style-type: none"> • 208–230V, 50–60 Hz Hz, 2A • If the Torsion power supply is integrated into the Axial power supply, rather than separate, the system usage rises to 12A
Hot/Cold Chamber	<ul style="list-style-type: none"> • 230V, 50–60 Hz Hz, 9.6A/2200W, 1 ph
Saline Bath	<ul style="list-style-type: none"> • 207–230V, 50–60 Hz Hz, 3A
Multiple Sample Fatigue Saline Chamber	<ul style="list-style-type: none"> • 207–230V, 50–60 Hz Hz, 2.7A
Thermoelectric Cooled Grips	<ul style="list-style-type: none"> • 100–240V, 50–60 Hz Hz, 2.85–1.19A • Chiller: 115–230V, 50–60 Hz, 4.5–2.3A
Power cords provided	<ul style="list-style-type: none"> • 6-20P plug for 230V systems • International: Line power cord provided is based on country



6-20P



Use power cords with plugs appropriate for your circuit.



Supply voltages lower than indicated may result in a degradation of performance.



Ensure that the mains assigned do not also supply power to noise generating equipment nearby, such as motors, welders, transformers, etc.



An independent heavy GROUND wire must be provided through the power hookup. Improper grounding may cause severe damage for which the supplier will not accept responsibility. All power strips must be fully grounded and carry the ground through to the sockets into which the computer is plugged.

Utility Requirements



GAS

Item	Requirement																																																																					
Purge gas	Air																																																																					
Pressure	413–515 kPa (60–80 psig)																																																																					
Flow	8.5 L/min (0.3 CFM)																																																																					
Connections	¼-inch push-to-connect or ¼-inch NPT female to filter/dryer/regulator with ¼-inch push-to-connect removed ¹																																																																					
Conditions	<table border="1"> <thead> <tr> <th rowspan="2">ISO8573-1:2010 Class</th> <th colspan="3">Solid Particulate</th> <th rowspan="2">Concentration mg/m³</th> <th colspan="2">Water</th> <th rowspan="2">Oil Total oil (aerosol, liquid, & vapor) ppm (mg/m³)</th> </tr> <tr> <th colspan="3">Maximum number of particles per m³</th> <th>Vapor</th> <th>Liquid</th> </tr> <tr> <td></td> <td>0.1 – 0.5 micron</td> <td>0.5 – 1 micron</td> <td>1 – 5 micron</td> <td></td> <td>Pressure dewpoint</td> <td>g/m³</td> <td></td> </tr> </thead> <tbody> <tr> <td>1</td> <td>≤ 20,000</td> <td>≤ 400</td> <td>≤ 10</td> <td>-</td> <td>≤ -94°F (-70°C)</td> <td>-</td> <td>0.008 (0.01)</td> </tr> <tr> <td>2</td> <td>≤ 400,000</td> <td>≤ 6,000</td> <td>≤ 100</td> <td>-</td> <td>≤ -40°F (-40°C)</td> <td>-</td> <td>0.08 (0.1)</td> </tr> <tr> <td>Recommended → 3</td> <td>-</td> <td>≤ 90,000</td> <td>≤ 1,000</td> <td>-</td> <td>≤ -4°F (-20°C)</td> <td>-</td> <td>0.83 (1)</td> </tr> <tr> <td>Minimum Specs → 4</td> <td>-</td> <td>-</td> <td>≤ 10,000</td> <td>-</td> <td>≤ +37°F (3°C)</td> <td>-</td> <td>4.2 (5)</td> </tr> <tr> <td>5</td> <td>-</td> <td>-</td> <td>≤ 100,000</td> <td>-</td> <td>≤ +45°F (7°C)</td> <td>-</td> <td>-</td> </tr> <tr> <td>6</td> <td>-</td> <td>-</td> <td>-</td> <td>≤ 5</td> <td>≤ +50°F (10°C)</td> <td>-</td> <td>-</td> </tr> </tbody> </table>	ISO8573-1:2010 Class	Solid Particulate			Concentration mg/m ³	Water		Oil Total oil (aerosol, liquid, & vapor) ppm (mg/m ³)	Maximum number of particles per m ³			Vapor	Liquid		0.1 – 0.5 micron	0.5 – 1 micron	1 – 5 micron		Pressure dewpoint	g/m ³		1	≤ 20,000	≤ 400	≤ 10	-	≤ -94°F (-70°C)	-	0.008 (0.01)	2	≤ 400,000	≤ 6,000	≤ 100	-	≤ -40°F (-40°C)	-	0.08 (0.1)	Recommended → 3	-	≤ 90,000	≤ 1,000	-	≤ -4°F (-20°C)	-	0.83 (1)	Minimum Specs → 4	-	-	≤ 10,000	-	≤ +37°F (3°C)	-	4.2 (5)	5	-	-	≤ 100,000	-	≤ +45°F (7°C)	-	-	6	-	-	-	≤ 5	≤ +50°F (10°C)	-	-
	ISO8573-1:2010 Class		Solid Particulate				Concentration mg/m ³	Water		Oil Total oil (aerosol, liquid, & vapor) ppm (mg/m ³)																																																												
		Maximum number of particles per m ³			Vapor	Liquid																																																																
		0.1 – 0.5 micron	0.5 – 1 micron	1 – 5 micron		Pressure dewpoint	g/m ³																																																															
	1	≤ 20,000	≤ 400	≤ 10	-	≤ -94°F (-70°C)	-	0.008 (0.01)																																																														
	2	≤ 400,000	≤ 6,000	≤ 100	-	≤ -40°F (-40°C)	-	0.08 (0.1)																																																														
	Recommended → 3	-	≤ 90,000	≤ 1,000	-	≤ -4°F (-20°C)	-	0.83 (1)																																																														
Minimum Specs → 4	-	-	≤ 10,000	-	≤ +37°F (3°C)	-	4.2 (5)																																																															
5	-	-	≤ 100,000	-	≤ +45°F (7°C)	-	-																																																															
6	-	-	-	≤ 5	≤ +50°F (10°C)	-	-																																																															

¹ Interconnection between the system and the supplied filter/regulator equipment is via ¼-inch push-to-connect fittings; use supplied tubing.



HOT/COLD CHAMBER

Item	Requirement
Gas	Liquid nitrogen
Pressure	152–345 kPa (22–50 psig)
Connections	½ -inch SAE, 45 degree flare fitting

TA Instruments Offices

For information on our latest products, contact information, and more, see our website at:
<http://www.tainstruments.com>.

To find your local TA Instruments office and contact information, visit
<http://www.tainstruments.com/contact/ta-directory/>

TA Instruments – Waters LLC
Corporate Headquarters
159 Lukens Drive
New Castle, DE 19720
USA

Telephone: 302-427-4000
Fax: 302-427-4001
Email: info@tainstruments.com